

FREQUENTLY ASKED QUESTIONS

Do I need to clean the surface or equipment before using Sani-Cloth® 70 to disinfect?

Yes, it is important that you clean the surface before disinfecting if it is visibly dirty. Sani-Cloth® Detergent is great for cleaning prior to disinfection with Sani-Cloth® 70.

What is the composition of Sani-Cloth® 70?

Active Ingredients: Isopropyl Alcohol 70% v/v.

Is Sani-Cloth® 70 effective against SARS-CoV-2?

Yes. SARS-CoV-2 is the name given to a strain of Coronavirus which is in the family of enveloped viruses. Vaccinia Virus is also within the family of enveloped viruses. This is perceived as more difficult to kill, therefore, any product with proven efficacy against vaccinia virus EN 14476 is by default effective against SARS-CoV-2 including all other mutations of the virus and any other enveloped viruses. Sani-Cloth® 70 contact time (for Clean and Dirty conditions) for Vaccinia Virus - Strain Elstree is 15sec.

What is the contact time and what happens if the surface dries before the states contact time in Sani-Cloth® 70 product labels?

The contact times listed in the efficacy table on the product label is the amount of time that surfaces should remain visibly wet to inactivate the listed microorganisms. The table presented shows the contact time and the test condition. Ensure the surface remains visibly wet for the required contact time.

In certain geographies and in settings where temperature, relative humidity and air changes may vary, it is possible that users may need to use an additional wipe to achieve the required contact time.

What surfaces and equipment are compatible with Sani-Cloth® 70?

Sani-Cloth® 70 wipes are designed to be compatible with medical devices: hard non-porous surfaces. Follow medical device manufacturer cleaning and disinfecting guidelines. If no manufacturers instructions for use are available, test in an inconspicuous area before use.

Are Sani-Cloth® 70 wipes a hazard to my patients and staff?

When used according to the label instructions Sani-Cloth® 70 should not present a safety hazard. Sani-Cloth 70® formulation does not contain any substances that are assessed to be PBT (persistent, bio accumulative and toxic) or vPvB (very persistent and very bio accumulative).

What PPE should I wear when using Sani-Cloth® 70 wipes?

We recommend wearing gloves when using Sani-Cloth® 70 wipes. Some other items of PPE will be determined by local guidelines and infection risks assessment.

The use of PPE will depend on the risk assessment and the microorganisms that you may come into contact with when cleaning and disinfecting, e.g. If you are cleaning and disinfecting equipment after patient use in an isolation room, your PPE requirements may be different to those required when cleaning and disinfecting a general use ward area. PPE is required to protect Healthcare Professionals from coming into contact with microorganisms when cleaning and disinfecting. Please follow local Infection Control protocols at all times on the use of PPE.

What is the expiration date of Super Sani-Cloth® 70 wipes?

Sani-Cloth® 70 wipes have a shelf life of 36 months from the date of manufacture. The expiration date is stamped directly onto the product packaging.

Once open how long do I have to use the product?

Sani-Cloth® 70 wipes have a Period After Opening of 3 months. Ensure the lid on the canister is closed when not in use, do not use any wipes which have become dry. Do not re-use any wipes.

How should Sani-Cloth® 70 be stored?

Sani-Cloth® 70 wipes should be stored in an upright position away from direct sunlight, in a dry cool, well ventilated area. Temperature 0°C - 25°C.

How should I dispose of Sani-Cloth® 70 ?

Follow your local healthcare setting waste disposal protocols.

Are Sani-Cloth® 70 wipes flushable or maceratable?

No, Sani-Cloth® 70 wipes are not flushable nor maceratable. The DO NOT FLUSH symbol is visible on the front of the packaging.

Can I use Sani-Cloth® 70 wipes if the pack is damaged or the lid is broken?

Sani-Cloth® 70 wipes should not be used if the pack is damaged or partially open in any way.

If I require some support or training on using this product, how do I get in touch with the company?

Please refer to the contact details below.



For more information, write to us: contact@pdi-emea.com

T: +44 (0) 8081 697 945

E: contact@pdi-emea.com

W: pdi-intl.com

 **PDI International**