



Sani-Cloth® AF: Product Guidance

Two in one cleaning & disinfectant alcohol-free wipes

For non-invasive medical devices and non-porous hard surfaces. It can be used to pre-clean or decontaminate medical devices prior to sterilisation.
For professional use only.

Preparation



Dispense



- + Dispense the wipe at a 45° angle, pulling away from you.
- + Close the lid when not in use.

Clean



- + If the surface is visibly soiled, clean the surface following an S-shaped motion from top to bottom.
- + Use additional wipes as needed.
- + Use a fresh wipe for each new surface.

Disinfect



- + Using a fresh wipe, disinfect the surface in an S-shaped motion from top to bottom, avoiding cross-contamination.
- + Use a fresh wipe for each new surface.

Contact Time



- + Ensure the surface remains visibly wet for 2 minutes for complete disinfection.
- + Use an additional wipe if necessary.
- + Allow to air dry.

Disposal



- + Dispose of wipe after a single use.
- + Packaging is not recyclable.

Other:
Printed version available on request



Sani-Cloth® AF: Product Guidance

Two in one cleaning & disinfectant alcohol-free wipes

Other information:

- + Store upright in a cool, dry place.
- + Avoid contact with eyes, skin & clothes.
- + Refer to the expiry date on the packaging.
- + Approx.1 wipe per 1m².
- + Follow local healthcare policy for personal protective equipment, hand hygiene & waste disposal guidelines.

Contra-indications:

- + Do not use for terminal disinfection of invasive medical devices.
- + Any serious incident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State.
- + Compatibility with medical devices should be checked before use



Product Efficacy

Sani-Cloth® AF	Test Method	Organism	Contact time	
Bacteria	EN 16615	Pseudomonas aeruginosa	60 sec (c)	2 min (d)
	EN 16615	Staphylococcus aureus	60 sec (c)	2 min (d)
	EN 16615	Enterococcus hirae	60 sec (c)	2 min (d)
	EN 13727	Pseudomonas aeruginosa	30 sec (c)	
	EN 13727	Staphylococcus aureus	30 sec (c)	
	EN 13727	Enterococcus hirae	30 sec (c)	
Yeast	EN 16615	Candida albicans	60 sec (c)	2 min (d)
	EN 13624	Candida albicans	30 sec (c)	
	EN 13624	Candida auris	30 sec (c)	
Enveloped virus	EN 14476	Vaccinia Virus	30 sec (c)	60 s (d)
	EN 14476	Corona,HIV, Hep & C, SARS, Flu	30 sec (c)	2 min (d)
Non-enveloped virus	EN 14476	Norovirus	30 sec (c)	60 sec (d)
TB	EN 14348	Mycobacterium terrae	2 min (c)	

**Clean Conditions (c) Dirty Conditions (d)*



For more information, write to us: contact@pdi-emea.com

Pywell Road, Willowbrook East Industrial Estate, Corby NN17 5XJ, UK

+ 44 (0) 8081 697 945

contact@pdi-emea.com

pdi-intl.com



PDI International