


Super Sani-Cloth®: Instructions for Use

Two in one cleaning & disinfectant wipes

Highly compatible disinfection wipe for non-invasive medical devices and non-porous hard surfaces. It can be used to pre-clean or decontaminate medical devices prior to sterilisation. For professional use only.



Preparation




Dispense



- + Dispense the wipe at a 45° angle, pulling away from you.
- + Close the lid when not in use.

Clean



- + If the surface is visibly soiled, clean the surface following an S-shaped motion from top to bottom.
- + Use additional wipes as needed.
- + Use a fresh wipe for each new surface.

Disinfect



- + Using a fresh wipe, disinfect the surface in an S-shaped motion from top to bottom, avoiding cross-contamination.
- + Use a fresh wipe for each new surface.



Contact Time



- + Ensure the surface remains visibly wet for 1 minute for complete disinfection.
- + Use an additional wipe if necessary.
- + Allow to air dry.

Disposal



- + Dispose of wipe after a single use.
- + Packaging is recyclable



Canister

Lid

Other:
Printed version available on request

Super Sani-Cloth®: Instructions for Use

Two in one cleaning & disinfectant wipes

Caution and other information:

- + IF IN EYES: Rinse cautiously for several minutes. If irritation persists, get medical advice
- + Highly flammable liquid and vapour
- + Keep product closed and stored upright in a cool dry place away from heat sources
- + Vapour may cause drowsiness or dizziness
- + Avoid release to the environment. Harmful to aquatic life with long-lasting effects
- + Refer to the expiry date on the packaging
- + Approx.1 wipe per 1m².
- + Dispose of contents and container in accordance with local regulations



DANGER Contains Propan-2-ol

Contra-indications:

- + Do not use for terminal disinfection of invasive medical devices.
- + Any serious incident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State.

Formulation

62%v/v Isopropyl alcohol (IPA) (55%w/w)

0.5% Dual Quaternary Ammonium Compounds (Quats)

Product Efficacy

Super Sani-Cloth®	Test Method	Organism	Contact time	
Bacteria	EN 16615	Pseudomonas aeruginosa	60 sec (c)	60 sec (d)
	EN 16615	Staphylococcus aureus	60 sec (c)	60 sec (d)
	EN 16615	Enterococcus hirae	60 sec (c)	60 sec (d)
	EN 13727	Pseudomonas aeruginosa	30 sec (c)	60 sec (d)
	EN 13727	Staphylococcus aureus	30 sec (c)	60 sec (d)
	EN 13727	Enterococcus hirae	30 sec (c)	60 sec (d)
Yeast	EN 16615	Candida albicans	60 sec (c)	60 sec (d)
	EN 13624	Candida albicans	30 sec (c)	
	EN 13624	Candida auris	30 sec (c)	
Enveloped virus	EN 14476	Vaccinia Virus	60 sec (c)	60 s (d)
	EN 14476	Corona, HIV, Hep & C, SARS, Flu	30 sec (c)	2 min (d)
Non-enveloped virus	EN 14476	Norovirus	30 sec (c)	60 sec (d)
TB	EN 14348	Mycobacterium terrae	60 sec (c)	

*Clean Conditions (c) Dirty Conditions (d)



For more information, write to us: contact@pdi-emea.com



PDI Ltd. Pywell Road, Willowbrook East Industrial Estate, Corby NN17 5XJ, UK

T: + 44 (0) 8081 697 945

E: contact@pdi-emea.com

W: pdi-intl.com



PDI International