


Super Sani-Cloth®: Product Guidance

Two in one cleaning & disinfectant wipes

Highly compatible disinfection wipe for non-invasive medical devices and non-porous hard surfaces. It can be used to pre-clean or decontaminate medical devices prior to sterilisation. For professional use only.

Preparation

Dispense

- + Dispense the wipe at a 45° angle, pulling away from you.
- + Close the lid when not in use.

Clean

- + If the surface is visibly soiled, clean the surface following an S-shaped motion from top to bottom.
- + Use additional wipes as needed.
- + Use a fresh wipe for each new surface.

Disinfect

- + Using a fresh wipe, disinfect the surface in an S-shaped motion from top to bottom, avoiding cross-contamination.
- + Use a fresh wipe for each new surface.

Contact Time

- + Ensure the surface remains visibly wet for 1 minute for complete disinfection.
- + Use an additional wipe if necessary.
- + Allow to air dry.

Disposal

- + Dispose of wipe after a single use.
- + Packaging is recyclable

Super Sani-Cloth®: Product Guidance

Two in one cleaning & disinfectant wipes

Caution and other information:

- + IF IN EYES: Rinse cautiously for several minutes. If irritation persists, get medical advice
- + Highly flammable liquid and vapour
- + Keep product closed and stored upright in a cool dry place away from heat sources
- + Vapour may cause drowsiness or dizziness
- + Avoid release to the environment. Harmful to aquatic life with long-lasting effects
- + Refer to the expiry date on the packaging
- + Approx. 1 wipe per 1m².
- + Dispose of contents and container in accordance with local regulations

DANGER Contains Propan-2-ol

Contra-indications:

- + Do not use for terminal disinfection of invasive medical devices.
- + Any serious incident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State.

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For more information, write to us: contact@pdi-emea.com

PDI Ltd. Pywell Road, Willowbrook East Industrial Estate, Corby NN17 5XJ, UK

T: + 44 (0) 8081 697 945

E: contact@pdi-emea.com

W: pdi-intl.com

PDI International