

Sani-Cloth® Detergent wipes are impregnated with a low-level detergent for the cleaning & damp dusting of all general surfaces and equipment, including beds and mattresses. They contain a clean, fresh fragrance with odour-neutralising properties, utilising a low-smear, biodegradable formulation.



Recommended Uses:

- + Routine cleaning procedures
- + Prior to Disinfecting in a two-step cleaning & disinfecting process

+ Ready to Use

Pre-dosed wipes requiring no mixing or measuring, with low smear formulation.

+ Free From

Alcohol free, bleach free and PHMB free, causing no skin irritation or sensitisation.

+ Supporting Resources

Point of Use accessories, "I Am Clean" tape & sticky notes, IFU posters.

+ Improved Efficiencies

Easy and single wipe dispensing saves valuable time.

+ Maximum Flexibility

Suitable for any healthcare environment available in soft pack, canisters and buckets.

+ Supports Carbon Net Zero Goals

Manufactured in the UK, resulting in supply chain reliability and carbon footprint reduction.

Product Details

Format	Wipe Size	Wipe Count	QTY per Case	Classification	Product code
225 ct Bucket	270mmx 267mm	225	1	MD Class IIa	EDT1210C
100 ct Canister	240mm X 200mm	100	6	MD Class IIa	EDT1230C
150 ct Canister	300mm X 200mm	150	6	MD Class IIa	EDT1240
100ct flow wrap	300mm X 200mm	100	12	MD Class IIa	EDT1250C



For more information, write to us: contact@pdi-emea.com

T: + 44 (0) 8081 697 945 E: contact@pdi-emea.com W: pdi-intl.com PDI International