

Easy Screen

CLEANING WIPE



Easy Screen® cleaning wipes contain 70% IPA, ideal for cleaning dirt, grime, and fingerprints on standard touch-screen healthcare equipment and more.

+ Ready to Use

Pre-dosed wipes requiring no mixing or measuring, with no smear formulation.

+ Optimal performance

Features the cleaning power of 70% Isopropyl alcohol (IPA).

+ Enhanced Cleaning

Fast drying, non-residue, anti-fogging, effectively removing streaks, smears and fingerprints.

+ Multi Surface Compatibility

Ideal for use on touch screens, tablets, electronic medical equipment, and procedural lights with substantial supportive compatibility data available.

+ Supporting Resources

Point of Use accessories, "I Am Clean" tape & sticky notes, IFU posters.

+ Supports Carbon Net Zero Goals

Manufactured in the UK, resulting in supply chain reliability and carbon footprint reduction.




Product Details

Format	Wipe Size	Wipe Count	QTY per Case	Classification	Product Code	NHS Code
100ct canister	200mm x 178mm	100	6	Consumable	EES1060	VJT595

Compatible with:

- + Polycarbonate
- + Corning® Gorilla® Glass
- + Sapphire glass
- + Aluminium silicate
- + Acrylic
- + Glass
- + Etched glass
- + Stainless steel
- + Plus more...



Before cleaning



After cleaning with Easy Screen



For more information, write to us: contact@pdi-emea.com

T: + 44 (0) 8081 697 945 E: contact@pdi-emea.com W: pdi-intl.com

 PDI International