

Easy Screen

CLEANING WIPE

WIPE



Easy Screen cleaning wipes contain 70% IPA, ideal for cleaning dirt, grime, fingerprints on common touch-screen healthcare equipment and more.

- + **Optimal performance:** features the cleaning power of 70% Isopropyl alcohol (IPA).
- + **Ready to use:** pre-dosed wipes requiring no mixing or measuring, with no smear formulation.
- + **Enhanced cleaning:** fast drying, non-residue, anti-fogging, effectively removing streaks, smears and fingerprints.
- + **Multi-surface compatibility:** ideal for use on touch screens, tablets, electronic medical equipment, and procedural lights with substantial supportive compatibility data available.



+ PDI
BE THE DIFFERENCE®



Touch-screens




Tablets, electronic medical record equipment



OR and procedural lights within a housing unit



Before cleaning



After cleaning with Easy Screen

- Fast-drying
- No residue
- Anti-fogging
- Anti-streaking

Compatible with

- Polycarbonate
- Corning® Gorilla® Glass
- Sapphire glass
- Aluminium silicate
- Acrylic
- Glass
- Etched glass
- Stainless steel

And more...

Point of Use Accessories



Bracket for canisters

Bracket: A001BLK

Canister brackets affix securely to the wall, holding the canister in place at its location.



Pole Adapter

Product Code: A002POL

Pole brackets make it possible to affix the canister brackets on to portable equipment such as BP stands.

Range Details

	REORDER NO	WIPE COUNT	WIPE SIZE	QTY PER CASE	CLASSIFICATION
100ct Large Canister	EES1060	100	200mm x 178mm	6	Consumable

Other formulations available:



Find out more about how our Layered Approach to Infection Prevention can help you.

PDI EMEA, Pywell Road, Willowbrook Est. Industrial Estate, Corby - Northants, NN17 5XJ, UK

+44 (0) 8081 697 945 contact@pdi-emea.com wearepdi-intl.com PDI International PDI_Intl

