

# **Prevantics®**

## **2% CHG Applicators**



**2% CHG  
APPLICATOR**

The single-use applicator is a ready to use sterile device containing 2% (w/w) Chlorhexidine gluconate and 70% (v/v) Isopropyl alcohol requiring no activation.

In the UK, these products are not licensed medicinal products, and they are not intended for pre surgical skin antiseptics. These products are PT1 Skin Biocide products, intended for skin disinfection only.



### **Ideal to use for:**

- + Skin Disinfection using 2% Chlorhexidine gluconate (CHG) and 70% Isopropyl alcohol (IPA).
- + Larger areas of skin that require disinfection.

## **Features & Benefits**

### **+ Proven Efficacy**

Based on its active ingredients, a broad spectrum of efficacy is guaranteed against bacterial, yeasticidal & virucidal activity.

### **+ Convenient to Use**

Ready to use, requiring no activation. Simply press the plunger to release the liquid onto the sponge allowing for effective skin disinfection.

### **+ Sterile**

The liquid has been filtered and then the applicator sterilised using Ethylene Oxide (ETO).

### **+ Time Saving**

With an integrated delivery system, and no sharp glass to break, the applicator avoids unwanted spillages and accidents.

### **+ Excellent User Experience**

The all-round foam sponge allows for gentle friction application improving user comfort.

### **+ Made in Europe**

We aim to provide peace of mind to our customers by guaranteeing compliance with rigorous regulations and standards, quality, uninterrupted supply, and reduced CO2 footprint.



## 2% CHG Applicators

### Recommended applicator based on general skin coverage



APPLICATOR SIZE	RECOMMENDED SKIN COVERAGE
26ml	50cm x 50cm
10.5ml	30cm x 30cm
3ml	15cm x 15cm

### Range Details

Format	Product Size	QTY per Case	Classification	Product code	NHS code	Route
26ml Clear Applicator	26ml	10 x 6 (60)	Biocide	IAP4140	DEC85017	NHS Supply Chain
26ml Tinted Applicator	26ml Tinted	10 x 6 (60)	Biocide	IAP4230	DMK85007	NHS Supply Chain
10.5ml Clear Applicator	10.5ml	10 x 6 (60)	Biocide	IAP4130	DMK070	NHS Supply Chain
10.5ml Tinted Applicator	10.5ml Tinted	10 x 6 (60)	Biocide	IAP4220	DMK85009	NHS Supply Chain
3ml Clear Applicator	3ml	10 x 6 (60)	Biocide	IAP4110	DMK85006	NHS Supply Chain
3ml Tinted Applicator	3ml Tinted	10 x 6 (60)	Biocide	IAP4200	DMK85010	NHS Supply Chain

*\*Tinted version available for training purposes*

#### Also available:

+ CHG swabsticks + CHG solutions + CHG full-body washcloths



**For more information, write to us: [contact@pdi-emea.com](mailto:contact@pdi-emea.com)**

T: + 44 (0) 8081 697 945

E: [contact@pdi-emea.com](mailto:contact@pdi-emea.com)

W: [pdi-intl.com](http://pdi-intl.com)

 PDI International