

# Wipe Selector Guide

## 2 in 1 Cleaning & Disinfecting



At PDI International, we understand there are different variables and sensitivities in every healthcare setting; **we have a solution to address every situation.** We can support you and your facility in recommending which wipe should be used to ensure patient and healthcare worker safety as well as equipment protection.

### Cleaning & Disinfection

- ✚ High touch surfaces in A&E, wards, waiting areas and operating theatres
- ✚ Mattresses, bedrails, sideboards, doorhandles
- ✚ Common equipment and non-invasive medical devices
- ✚ Busy Areas with high traffic; waiting rooms & porters chairs

#### Recommended Products:



#### Super Sani-Cloth®

Ideal for fast paced environments when quick drying is necessary



#### Sani-Cloth® AF

Ideal for surfaces where alcohol is not recommended

### Outbreaks & Isolation Situations

- ✚ Confirmed or suspected cases of Clostridium difficile or other known infections
- ✚ All surfaces including commodes and bathrooms

#### Recommended Product:



#### Sani-Cloth® Chlor

### Equipment Compatibility

- ✚ Fast drying with broad microorganism efficacy and compatibility
- ✚ Compatible with more than 600 manufacturer's IFUs including patient monitors, ventilators, keyboards, PC's, BP machines, hand held devices, and touch screens\*

#### Recommended Product:



#### Super Sani-Cloth®

Ideal for fast paced environments when quick drying is necessary



### Sani-Cloth® Chlor

#### Sporicidal

2-in-1 Cleaning & Disinfectant

- ✚ Complete disinfection in 2 minutes
- ✚ Product efficacy:
  - ✓ Sporicidal
  - ✓ Tuberculocidal
  - ✓ Virucidal
  - ✓ Yeasticidal
  - ✓ Bactericidal



Find Out More



### Super Sani-Cloth®

#### Alcohol & Quat Formulation

2-in-1 Cleaning & Disinfectant

- ✚ Complete disinfection in 1 minute
- ✚ Compatible by Design
- ✚ Product efficacy:
  - ✓ Tuberculocidal
  - ✓ Virucidal
  - ✓ Yeasticidal



Find Out More



### Sani-Cloth® AF

#### Alcohol Free

2-in-1 Cleaning & Disinfectant

- ✚ Complete disinfection in 2 minutes
- ✚ For surfaces where alcohol is not recommended
- ✚ Product efficacy:
  - ✓ Tuberculocidal
  - ✓ Virucidal
  - ✓ Yeasticidal
  - ✓ Bactericidal



Find Out More

\*Check with Manufacturer's IFUs



For details of our full range of surface cleaning & disinfection solutions contact us.



contact@pdi-emea.com



+44 (0) 8081 697 945



wearepdi-intl.com



PDI International



@pdi\_intl