



Prevantics 2% CHG Swabstick: Instructions For Use

Please follow local policy for personal protective equipment, hand hygiene & waste disposal guidelines.



Step 1: Safely open the packaging without any external contact with the product.



Step 2: Safely remove the swabstick from the packaging without any external contact.



Step 3: Apply the the product on the area to be disinfected using vigorous back-and-forth, up-and-down, repeated strokes.



Step 4: Turning the swabstick over half way through, for at least 30 seconds, allow the skin to air dry.

SWABSTICK SIZE	RECOMMENDED SKIN COVERAGE
1.6ml	10cm x 10cm



For more information, write to us: contact@pdi-emea.com

T: + 44 (0) 8081 697 945

E: contact@pdi-emea.com

W: wearepdi-intl.com

PDI International