

## 2% CHG Device Swab

## Device SWAB

Quick and easy to use, its innovative design and trusted formulation of Chlorhexidine gluconate (CHG) at 2% and Isopropyl alcohol (IPA) at 70% help standardise practice when disinfecting needle-free injection sites.

### Ideal to use for:

+ Medical devices such as external parts of needleless connectors.



### + Proven efficacy

A broad spectrum of efficacy guaranteed against bacterial, yeasticidal & virucidal activity.

### + Moisture retention

Airtight poly-lined foil sachet to minimise the risk of the swab drying.

### + Convenient to use

Individually wrapped single-use swabs with easy tearing.

## Efficacy Table

2% CHG SWAB	TEST METHOD	ORGANISM	CONTACT TIME*	
Bacteria	EN 16615	Pseudomonas aeruginosa	60 sec (c)	
	EN 16615	Staphylococcus aureus	60 sec (c)	
	EN 16615	Enterococcus hirae	60 sec (c)	
	EN 13727	Pseudomonas aeruginosa	30 sec (c)	
	EN 13727	Staphylococcus aureus	30 sec (c)	
	EN 13727	Enterococcus hirae	30 sec (c)	
Yeasts	EN 16615	Candida albicans	60 sec (c)	
	EN 13624	Candida albicans	30 sec (c)	
Enveloped Virus	EN 14476	Vaccinia Virus	15 sec (c)	60 sec (d)
	EN 14476	Corona, HIV, Hep B & C, SARS, Flu	15 sec (c)	60 sec (d)
TB	EN 14348	Mycobacterium terrae	60 sec (c)	

\*Clean Conditions (c) / Dirty Conditions (d)

## Product Details

PRODUCT SIZE	WIPE COUNT	QTY PER CASE	CLASSIFICATION	REORDER NO
150mm x 129mm	1	100 x 12 (1200)	MD Class IIa	ISW4040

2% Chlorhexidine gluconate (CHG) 70% Isopropyl alcohol (IPA) formulation is recommended in evidence-based infection guidelines published by health organisations such as NICE<sup>1</sup>, epic3<sup>2</sup>, and CDC<sup>3</sup>



For more information, write to us: [contact@pdi-emea.com](mailto:contact@pdi-emea.com)

T: + 44 (0) 8081 697 945 E: [contact@pdi-emea.com](mailto:contact@pdi-emea.com) W: [wearepdi-intl.com](http://wearepdi-intl.com) PDI International