

# Product Portfolio







## Contents

<b>Welcome to PDI International</b> .....	<b>3</b>
<b>Delivering Meaningful Impact</b> .....	<b>4</b>
<b>Why Partner with PDI</b> .....	<b>5</b>
<b>Infection Risks Are Everywhere, So Are We</b> .....	<b>6</b>
<b>Sustainable Innovation</b> .....	<b>7</b>
<b>Layered Approach to Infection Prevention</b> .....	<b>8-9</b>
<b>Product Range</b> .....	<b>10-11</b>
<b>Cleaning &amp; Disinfection</b> .....	<b>12-17</b>
<b>Automated Room Disinfection</b> .....	<b>18-19</b>
<b>Skin Disinfection</b> .....	<b>20-25</b>
<b>Patient Cleansing</b> .....	<b>26-27</b>

## Welcome to PDI International

### About Us

PDI International is dedicated to reduce infections and promote health and wellness in healthcare, hospitality and our communities.

Driven by a commitment to research and development, sustainable innovation and industry-leading partnerships, our range of environmental cleaning and disinfection, skin decontamination and patient care solutions are trusted by healthcare professionals around the world.

### Dedicated to Be The Difference®

The solutions provided by PDI International are designed to respond to the specific needs and regulatory requirements in each unique environment where we work. It is only by responding to those specific needs we can truly make a difference, and impact lives within each of those markets.

### Our Mission

We are driven by a higher purpose to lead the fight against preventable infections and help save lives.

Because of what we do every day, someone's Mother, Brother, Father, Son, Sister, Aunt, Uncle, Child or Lover won't get sick, die or suffer.

### Our Vision

To be a global infection prevention innovative leader and trusted advisor to our customers internationally.

Together, we will make all the difference in the world.



# Delivering Meaningful Impact

- 1.** The "Julius" family of PDI invented the very first wet wipe known as the "moist towelette" which was later provided with every KFC meal!
- 2.** We have developed some of the most important infection prevention innovations, including the first alcohol prep pad and disinfectant disposable wipe.
- 3.** We are market leaders in Healthcare. 250,000,000 people are helped by PDI products every year globally.
- 4.** We have been manufacturing private label wet wipe solutions for global healthcare companies for more than 30 years.
- 5.** Every summer, we celebrate Be The Difference® day which represents a collective commitment to meaningfully engage, educate, and give back to the communities.



It's all part of our passion to **Be The Difference®**, every day.

# Why Partner with PDI International



## We Are Family Owned

We are a family-owned company built on a 70-year legacy that remains authentically committed to protecting people through evidence-based solutions proven to reduce HCAs. We also place immense importance on engaging, educating and giving back meaningfully to communities, both locally and internationally.



## A Layered Approach to Infection Prevention

Our layered approach to infection prevention puts people at the heart of preventing infections. We like to understand the environment and challenges before advising on the most suitable product selection which is accompanied by continual training, education and monitoring compliance to ensure that the strongest line of protection is in place.



## We are Your Trusted Advisor

We will always do our best to find a solution to overcome customer challenges, but similarly we will be honest if we can't. By demonstrating this level of transparency, we have earned the privilege of being a trusted advisor to our customers internationally.



## Made in Europe


Our products are manufactured in Europe, at our facilities in the UK and Italy. With European manufacturing, we aim to provide peace of mind to our customers by guaranteeing compliance with rigorous regulations and standards, quality, uninterrupted supply, and reduced CO2 footprint.



# Infection risks are everywhere. So are we.

Our solutions are designed to reduce preventable infections and help save lives across the healthcare continuum.

## Community-based care



## Acute care




## Post-acute care



# Sustainable Innovation

To achieve the highest standards for our products and operations, we drive sustainable innovation in product design and manufacturing across the life cycle of our products.





# Layered Approach to Infection Prevention

No single approach can fully eliminate the risk of healthcare-associated infections (HAIs). That's why healthcare facilities need multiple layers of defence to truly combat infection. Our layered, consultative approach is designed to help you go beyond just product selection to focus on your overall infection prevention strategy using the following steps:

## 1. Assessment

We aim to understand your clinical setting requirements and challenges before putting forward the most suitable solution which will support you in achieving your specific infection prevention objectives.

## 2. Selection

Based on the outcome of the assessment, we will advise you which products across our three ranges will best overcome your challenges and support your goals.

## 3. Training

Effective, on-going training with access to a comprehensive resource library to ensure staff compliance.

## 4. Support

We are committed to delivering continuous training and education so that staff feel confident and informed on the latest infection prevention and control practices at all times.





# Product Range

PDI International offers a broad range of evidence-based, market-leading Environment, Surface, Skin and Patient Care solutions, all designed to help reduce preventable infections, control healthcare costs, and, ultimately, help save lives. It is all part of our passion to Be The Difference®, every day.

## Environment & Surface Care



### Cleaning & Disinfection

Surface cleaning and disinfectant solutions used for preventative purposes and outbreak situations.



### Automated Room Disinfection

Highly efficacious and versatile system which delivers on speed, safety, compatibility and cost.

**Sani-Cloth®**

**HDM**  
HyperDRY Mist

**Prevantics®**

**Hygea**

## Skin & Patient Care



### Skin Disinfection

Comprehensive range of skin and hand disinfectant solutions clinically effective against a wide range of pathogens.



### Patient Cleansing

Personal hygiene, cleansing and continence care solutions designed to make patients feel clean and comfortable.



# Cleaning & Disinfection



## Wipe Selector Guide

We understand there are different variables and sensitivities in every healthcare setting. We can support you and your facility in recommending which wipe should be used to ensure patient and healthcare worker safety as well as equipment protection.



	Sani-Cloth® Chlor	Super Sani-Cloth®	Sani-Cloth® AF	Sani-Cloth® 70%	Sani-Cloth® Detergent	Easy Screen
<b>Performance</b>						
Cleaning	✓	✓	✓		✓	✓
Disinfecting	✓	✓	✓	✓		
<b>Use</b>						
Routine between patients		✓	✓	✓	✓	✓
Outbreak prevention & management	✓	✓	✓	✓		
C.diff Outbreak	✓					
<b>Formulation</b>						
Sodium Hypochlorite	✓					
Dual quat + alcohol		✓				
Quat (alcohol-free)			✓			
Alcohol				✓		✓
Detergent					✓	
<b>Efficacy</b>						
Sporicidal - C.diff	✓					
Tuberculocidal		✓	✓	✓		
Virucidal - Norovirus	✓	✓	✓			
Virucidal - Enveloped Viruses	✓	✓	✓	✓		
Yeasticidal	✓	✓	✓	✓		
Bactericidal	✓	✓	✓	✓		
<b>Contact time</b>						
Total contact time of listed pathogens	2 minutes	1 minute	2 minutes	1 minute	N/A	N/A



## Sani-Cloth® Chlor

SPORICIDAL CLEANING AND DISINFECTANT WIPE



Sani-Cloth® Chlor is a sporicidal two-in-one cleaning and disinfectant solution ideal for use in high risk, contaminated environments. Sporicidal, virucidal, bactericidal, and yeasticidal in 2 minutes.

	WIPE SIZE	WIPE COUNT	QTY PER CASE	CLASSIFICATION	REORDER NO
50ct Large Canister	200mm x 195mm	50	6	MD Class IIa	ECH1400

## Sani-Cloth® 70%

DISINFECTANT WIPE



Sani-Cloth® 70 , containing 70% Isopropanol alcohol, is a multi-purpose solution ideal for the disinfection of non-invasive medical devices. Tuberculocidal, virucidal, bactericidal, and yeasticidal in 2 minutes.

	WIPE SIZE	WIPE COUNT	QTY PER CASE	CLASSIFICATION	REORDER NO
200ct Large Canister	200mm x 200mm	200	6	MD Class IIa	E701020

## Super Sani-Cloth®

CLEANING AND DISINFECTANT WIPE



Super Sani-Cloth® is a two-in-one dual quaternary ammonium and alcohol-based solution, ideal for cleaning and disinfecting in fast-paced environments that require short contact times and broad microbial coverage. "Purple top" is compatible in more than 700 manufacturer's IFU's. Tuberculocidal, virucidal, bactericidal, and yeasticidal in 1 minute.

	WIPE SIZE	WIPE COUNT	QTY PER CASE	CLASSIFICATION	REORDER NO
100ct Large Canister	200mm x 195mm	100	6	MD Class IIa	ESC1620

## Sani-Cloth® Detergent

DISPOSABLE CLEANING WIPE



Sani-Cloth® Detergent is a multi surface solution with a refreshing, light lemon fragrance, ideal for routine cleaning and damp dusting.

	WIPE SIZE	WIPE COUNT	QTY PER CASE	CLASSIFICATION	REORDER NO
225ct Bucket	270mm x 267mm	225	1	Consumable	EDT1210
150ct Large Canister	300mm x 200mm	150	6	Consumable	EDT1240
100ct Soft Pack	300mm x 200mm	100	12	Consumable	EDT1250

## Sani-Cloth® AF

CLEANING AND DISINFECTANT WIPE



Sani-Cloth® AF is an alcohol free, quaternary ammonium formulation providing a two-in-one solution ideal for cleaning and disinfecting. Tuberculocidal, virucidal, bactericidal, and yeasticidal in 2 minutes.

	WIPE SIZE	WIPE COUNT	QTY PER CASE	CLASSIFICATION	REORDER NO
100ct Soft Pack	280mm x 200mm	100	12	MD Class IIa	EAF1100
225ct Bucket	270mm x 267mm	225	1	MD Class IIa	EAF1110
200ct Soft Pack	267mm x 210mm	200	6	MD Class IIa	EAF1130
40ct Soft Pack	200mm x 200mm	40	24	MD Class IIa	EAF1140

## Easy Screen® 70%

CLEANING WIPE




Easy Screen® cleaning wipes contain 70% Isopropyl alcohol (IPA), ideal for cleaning dirt, grime, fingerprints on common touch-screen healthcare equipment and more.

	WIPE SIZE	WIPE COUNT	QTY PER CASE	CLASSIFICATION	REORDER NO
100ct Large Canister	200mm x 178mm	100	6	Consumable	EES1060

## Prevantics® 2% Device


SWAB



Ready to use disposable disinfection wipes for medical devices containing Chlorhexidine gluconate (CHG) at 2% and Isopropyl alcohol (IPA) at 70%.

	PRODUCT SIZE	WIPE COUNT	QTY PER CASE	CLASSIFICATION	REORDER NO
Device Swab	150mm x 129mm	1	100 x 12 (1200)	MD Class IIa	ISW4040





## Point of Use Accessories

Our point of use accessories help put PDI products in the right place at the right time to help you do your job more efficiently and effectively.



### Wall Mount Case for Flow Wrap

Our innovative design holds the lid of the flow wrap securely in place while in use. We use antimicrobial properties within the plastic itself, that is proven to reduce the growth of bacteria by up to 99.99% (ISO 22196:2011).



	COLOUR	REORDER NO	NHS CODE	ROUTE
Wall Mount Bracket for 200ct Flow Wrap	Green	A003GRN	N/A	Direct / Distributor

### Canister Brackets



Our canister brackets affix securely to the wall, they hold the canister in place using a twist lock motion provided by the design of our canisters, meaning once in place, the canister is held securely allowing easy, one handed wipe dispensing.

	COLOUR	REORDER NO
Sani-Care Brackets	Grey with Logo	A001GRY
Sani-Care Brackets	Yellow	A001YEL
Sani-Care Brackets	Red	A001RED
Sani-Care Brackets	Green	A001GRN
Sani-Care Brackets	Orange	A001ORG
Sani-Care Brackets	Purple	A001PUR
Sani-Care Brackets	Black	A001BLK

### I Am Clean Sticky Note

Our recyclable I am Clean sticky notes have been designed with users in mind. They can be used to easily and effectively identify equipment or surfaces that have been cleaned and stick easily to a clean, dry surface. They are paper based – meaning they can be written on using a ball point pen.



	COLOUR	COUNT	REORDER NO	NHS CODE	ROUTE
I am Clean - Sticky Notes	White and green	250 notes/pack (case of 96 packs)	A005SNT	VJT822	NHS Supply Chain

### I Am Clean Sticky Tape

The I Am Clean tape can be dispensed to whatever length is needed to affix to surfaces or bridge across gaps to clearly identify and communicate that a surface or a piece of equipment has been cleaned. There is space on the tape to enter the date, time and signature of the last clean.



	COLOUR	COUNT	REORDER NO	NHS CODE	ROUTE
I am Clean - Tape	White and green	Case of 72 rolls per case of 96 packs)	A006TPE	VJT838	NHS Supply Chain



# Automated Room Disinfection



HyperDRYMist® automated room disinfection is a highly efficacious and versatile system which provides excellent medical device compatibility, increased safety with low environmental impact.



## Micro-Nebulizer

The portable and lightweight device converts the disinfectant solution into an ultra-fine dry mist solution (95% particles < 1µm), resulting in complete and homogeneous coverage, delivering exceptional material compatibility, safely. Electronic programming with delayed start option and treatment reporting software available for compliance.



	DEVICE SIZE	QTY	CLASSIFICATION	REORDER NO
99MB Micro-Nebulizer	25cm x 42cm x 50cm	1	Low Voltage Equipment MD Class I	99MB

## Disinfectant Solution

Proprietary formula containing a patented non-toxic formulation containing low level hydrogen peroxide (6.6%) with carefully selected co-formulants, boosting its biocidal activity and overall stability making it effective against spores, fungi, viruses, bacteria and mycobacteria.




	DEVICE SIZE	QTY	CLASSIFICATION	REORDER NO
99S Disinfection Solution	1000ml	6	MD Class IIa	99S

For more information, request a copy of our HyperDRYMist® brochure.



# Skin Disinfection



## Prevantics® 2% CHG APPLICATOR

This single-use product contains filtered Chlorhexidine gluconate (CHG) at 2% w/w and isopropyl alcohol (IPA) at 70% v/v within a sterile applicator, providing a wide spectrum of efficacy against bacteria, mycobacteria, fungi and viruses. Available as clear or tinted solutions\*.



	PRODUCT SIZE	QTY PER CASE	CLASSIFICATION	REORDER NO
26ml Clear Applicator	26ml	10 x 6 (60)	Biocide	IAP4140
26ml Tinted Applicator	26ml Tinted	10 x 6 (60)	Biocide	IAP4230
10.5ml Clear Applicator	10.5ml	10 x 6 (60)	Biocide	IAP4130
10.5ml Tinted Applicator	10.5ml Tinted	10 x 6 (60)	Biocide	IAP4220
3ml Clear Applicator	3ml	11 x 6 (66)	Biocide	IAP4110
3ml Tinted Applicator	3ml Tinted	11 x 6 (66)	Biocide	IAP4200

\* Tinted versions available for training purposes

## Prevantics® 2% CHG SWABSTICK

The ready to use, premoistened swabstick contains 2% w/w Chlorhexidine gluconate (CHG) and 70% v/v Isopropyl alcohol (IPA) and has a unique 2-sided applicator for interdigital areas and skin folds ensuring thorough skin disinfection.



	PRODUCT SIZE	QTY PER CASE	CLASSIFICATION	REORDER NO
1.6ml Swabstick	1.6 ml	50 x 10 (500)	Biocide	ISS4010

## Prevantics® 2% CHG SWAB

The ready to use, premoistened wipe contains 2% w/w Chlorhexidine gluconate (CHG) and 70% v/v Isopropyl alcohol (IPA) to help prevent the spread of infections and help improve patient outcomes.



	PRODUCT SIZE	QTY PER CASE	CLASSIFICATION	REORDER NO
Skin Swab	150mm x 129mm	200 x 8 (1600)	Biocide	ISW4030

## Prevantics® RED 2% CHG SOLUTION



Ready to use skin disinfection biocide solution containing Chlorhexidine gluconate (CHG) at 2% w/w and Isopropyl alcohol (IPA) at 70% v/v providing long-lasting microbial kill of up to 48 hours. With optional red dye for superior visibility.

	PRODUCT SIZE	QTY PER CASE	CLASSIFICATION	REORDER NO
Bottle 120ml solution with optional red dye	120ml	6 x 8 (48)	Biocide	IRD6060

## Prevantics® 2% CHG SOLUTION



Ready to use skin disinfection biocide solutions containing Chlorhexidine gluconate (CHG) at 2% w/w and Isopropyl alcohol (IPA) at 70% v/v providing long-lasting microbial kill of up to 48 hours with either screw cap lid or applicator nozzle for accurate application.

	PRODUCT SIZE	QTY PER CASE	CLASSIFICATION	REORDER NO
Clear solution with applicator nozzle	200ml	24	Biocide	ILQ6030
Clear solution with applicator nozzle	500ml	24	Biocide	ILQ6050
Clear solution with screw top	500ml	24	Biocide	ILQ6040

## Prevantics® 4% CHG SOLUTION



Our CHG skin antiseptic contains a 4% w/w Chlorhexidine (CHG) cleansing formulation. CHG has a broad range of efficacy and is proven to be effective at reducing bacterial counts on the skin at the time of use, and for a period of time afterwards of up to 24 hours.

	PRODUCT SIZE	QTY PER CASE	CLASSIFICATION	REORDER NO
Bottle 500ml CHG 4% Liquid	500ml	20	Biocide	ILQ6000
Bottle 250ml CHG 4% Liquid	250ml	34	Biocide	ILQ6020
Bottle 100ml CHG 4% Liquid	100ml	40	Biocide	ILQ6010

## Prevantics® 70% IPA PREP PAD



The ready to use, premoistened wipe contains 70% Isopropyl alcohol (IPA) and water. Swabs are for skin preparation prior to injection.

	PRODUCT SIZE	QTY PER CASE	CLASSIFICATION	REORDER NO
Alcohol Prep Pad	65mm x 30mm	100 x 40 (4000)	MD Class I	IPP4060

## Prevantics® 2% CHG WASHCLOTHS



Large antiseptic body washcloths containing 2% Chlorhexidine (CHG), lightly fragranced, rinse-free providing rapid bactericidal action and long lasting residual protection. These wipes can be used straight from the pack or briefly warmed in the microwave prior to use.

	PRODUCT SIZE	WIPE COUNT	QTY PER CASE	CLASSIFICATION	REORDER NO
Full Body Washcloths	300mm x 200mm	8	36	Biocide	IWC3030

## Prevantics® 4% CHG KIT



The shower bath kit is suitable for pre-admission use to reduce the micro-organism count on the skin. Can be used in the shower or bath. The kit contains 100ml of 4% CHG solution and 2 dry sponges. The packaging provides easy to follow instructions on how to use.

	PRODUCT SIZE	QTY PER CASE	CLASSIFICATION	REORDER NO
Shower bath kit	2 x dry sponge and 100ml bottle	42	Biocide	ISK3040

## Prevantics® CHG SPONGE



Single-use bathing sponges containing Chlorhexidine gluconate (CHG) providing an antiseptic solution for skin cleansing and disinfection with a broad spectrum of efficacy.

	PRODUCT SIZE	QTY PER CASE	CLASSIFICATION	REORDER NO
4% CHG Sponge Tinted	120mm x 75mm	80	Biocide	ISP3000
2% CHG Sponge Tinted	120mm x 75mm	80	Biocide	ISP3010
2% CHG Sponge Clear, No Rinse	120mm x 75mm	80	Biocide	ISP3020






**Prevantics®** **4% CHG**  
BRUSH-SPONGE



Our CHG Brush-Sponge contains an antiseptic cleansing formulation of 4% w/w Chlorhexidine enhanced with advanced lathering agents and emollients. The ergonomic design allows for intuitive use of the soft sponge on one side and the flexible bristles on the other for effective cleansing. Includes nail pick.

	QTY PER CASE	CLASSIFICATION	REORDER NO
4% CHG Brush-Sponge with nail pick	25 x 9 (225)	MD Class I	IBS5000

**Brush-Sponge Wall Dispenser**



Our Brush-Sponge dispenser fixes securely to the wall. The stainless steel bracket holds the carton of 25 Brush-Sponges. The design of the bracket ensures enough packaging is visible for easy identification and the carton dispenser opening allows for quick product access.

	REORDER NO
Brush-Sponge Wall Dispenser	A004STL

**Prevantics®** **7.5% PVP**  
BRUSH-SPONGE



Our Iodine Brush-Sponge contains an antiseptic cleansing formulation of 7.5% Povidone Iodine (0.75% free iodine). Povidone Iodine has a broad spectrum of antimicrobial activity. The ergonomic design allows for intuitive use of the soft sponge on one side and the flexible bristles on the other for effective cleansing. Includes nail pick.

	QTY PER CASE	CLASSIFICATION	REORDER NO
7.5% Iodine Brush-Sponge with nail pick	25 x 9 (225)	MD Class I	IBS5040

**Hygea Hand Wipe**



Hygea Hand Wipes have been developed for quick hand cleansing and disinfection when there is limited access to soap and water. Clinically proven to kill 99.999% of bacteria. Effective against non-enveloped viruses including Norovirus in 30 seconds. Providing the most effective and hygienic alternative to soap and water.

	WIPE SIZE	WIPE COUNT	CASE QTY	CLASSIFICATION	REORDER NO
Single Sachet	175mm x 150mm	1	800	Biocide	PHW5500E
12ct Soft pack	175mm x 150mm	12	42	Biocide	PHW5520E
24ct Soft pack	175mm x 150mm	24	24	Biocide	PHW5530E
200ct Canister	195mm x 135mm	200	6	Biocide	PHW5540E

**Prevantics®** **DRY D1**  
BRUSH-SPONGE



Our DRY Brush-Sponge is sterile. The ready to use Brush-Sponge is compatible with any kind of antiseptic soaps and mild soaps. The product contains no detergent or antiseptic. The ergonomic design allows for intuitive use of the soft sponge on one side and the flexible bristles on the other for effective cleansing. Includes nail pick.

	QTY PER CASE	CLASSIFICATION	REORDER NO
Dry, Sterile Brush-Sponge with nail pick	25 x 9 (225)	MD Class IIa	IBS5020





## Hygea

### Hygea Continance Care PERSONAL WASHCLOTHS



Hygea Continance Care Washcloths containing 4% dimethicone are an all-in-one gentle cleansing and moisturising solution that provide barrier protection against Incontinence-Associated Dermatitis (IAD).

	WIPE SIZE	WIPE COUNT	CASE QTY	CLASSIFICATION	REORDER NO
24ct Soft Pack - Fragranced	300mm x 200mm	24	12	Cosmetic	PCC7120
8ct Soft Pack - Fragranced	300mm x 200mm	8	36	Cosmetic	PCC7110

### Hygea Full Body WASHCLOTHS



A waterless bed bathing solution pre-impregnated with skin cleansers and moisturisers to gently bathe and hydrate the skin. These washcloths can be used straight from the pack or can be briefly warmed in the microwave prior to use.

	WIPE SIZE	WIPE COUNT	CASE QTY	CLASSIFICATION	REORDER NO
8ct Soft Pack - Fragranced	300mm x 200mm	8	36	Cosmetic	PWC7100

### Hygea Hands & Face PERSONAL WASHCLOTHS



Hygea Hands & Face are patient care washcloths providing optimum skin cleansing, that are kind to the skin, leaving it clean and fresh.

	WIPE SIZE	WIPE COUNT	CASE QTY	CLASSIFICATION	REORDER NO
80ct Soft Pack - Fragranced - Maceratable	195mm x 155mm	80	12	Cosmetic	PHF7020





PDI EMEA Ltd. Pywell Rd, Willowbrook East Industrial Estate, Corby NN17 5XJ UK



+44 (0) 8081 697 945



contact@pdi-emea.com



pdi-intl.com